

MICRO CENTRIFUGE PURISPIN 6

컴팩트한 디자인의
개인용 미니 원심분리기



STANDARD SPECIFICATION

Max. RPM	6,200RPM (speed adjustments unavailable)	
Max. RCF	1,800 xg	
Max. Capacity	8 x 0.2/0.5/1.5/2.0ml tubes, 4 x 0.2ml PCR tube strips	
Rotor	Angel Type	8 Holes Circular Fixed-Angel rotor for 0.5/1.5/2.0 Tubes
	Strip Type	32 Hole Strip Rotor for 0.2 PCR Strip or individual Tubes
Speed Control	Automatic/Constant speed	
Driving Motor	DC motor	
Permissible Temp. & Humidity	5°C~40°C, <80%	
Rotating Direction	Clockwise	
Adaper	0.2/0.5ml	
Power Consumption	8W	
Safety Device	Door Limit Switch	
Metrial	Lid	PC (Polycarbonate)
	Body	ABS
	Rotor	ABS
Deimension (W x D x H) & Net Weight	155 x 183 x 125mm, 1.34Kg	
Packing Size (W x D x H) & Gross Weight	190 x 250 x 215mm, 1.60Kg	
Power Supply	1 Phase, AX 120V, 60Hz or AC 230V, 50/60 Hz	

MAIN FEATURE



컴팩트한 구조



우수한 내구성



심플한 작동법



저소음 구현

ACCESSORIES



Fixed angle rotor



Fixed angle rotor
for PCR strip



Adaptors for 0.2ml
and 0.4/0.5ml; tubes

ORDER INFORMATION

ITEM	DESCRIPTIONS	ORDER CODE
PURISPIN 6 Set	PURISPIN 6, with rotor Fixed Angle Rotor , Fixed Angle Rotor for PCR strip, Adaptors for 0.2ml and 0.4/0.5ml tubes	P6 - SET