

CRYSTE

MULTI CENTRIFUGE VARISPIN 4A

자주 사용하는 어플리케이션을 위한
앵글로터 전용 소형 원심분리기



STANDARD SPECIFICATION



Max. RPM	4,000rpm
Max. RCF	2,075 xg
Max. Capacity	15ml x 8
Accel Decel time to max rpm	20sec / 20sec
Noise Level	≤60dB
Set Time	99 min 59 sec or continuous
Time Count	From start / From set-up RPM
Program	10 memory
Display	Blue LCD
RPM Range	400 to 4,000 RPM
Display Message on screen	RPM(RCF), Time, Program
Available Rotor Quantity	1
Power Supply	V4A-MB
	AC 230V ± 10%, 50-60Hz, 1Phase
Power Consumption	0,15KW
Dimension(W x D x H)	291 x 416 x 244mm(H:550mm with open lid)
Weight (without rotor)	14kg

MAIN FEATURE



컴팩트한 구조

동급 최소 수준의 바닥면적으로
여러 장비를 사용하는 복잡한 실험테이블 위에서,
보다 자유로운 공간배치가 가능합니다.



높은 시인성의 LCD 스크린 사용



SOFT CLOSING DOOR LOCK



저소음 구현

AVAILABLE ANGLE ROTORS

AL815D	F15A-3	F15A-5				
 Fixed Angle rotor $\angle 40^\circ$ 8 x 15ml Radius(mm) : 116						
						
Volume (ml)	3	5	10	15	15 Conical	
Number per rack / rotor	8	8	8	8	8	
Adaptor bore (mm)	13.5	13.5	18.0	18	18	
Max. height for use (mm)	120	142	131	131	131	
RPM	4,000	4,000	4,000	4,000	4,000	
Radius (mm)	97.4	114.3	116	116	110.2	
RCF	1,742	2,045	2,075	2,075	1,971	

ORDER INFORMATION

	ITEM	DESCRIPTIONS	ORDER CODE
Main Body	VARISPIN 4A	VARISPIN 4A, without rotor, 230V/50~60Hz	V4A-MB
Angle Rotor	AL815D	15ml x 8 fixed angle rotor with 15ml stainless steel sleeves (Max 4,000RPM)	V4A-AR1