

VORTEX MIXER
PURIMIX



액체 시료의 간단한 믹싱을 위한
초소형 볼텍스 믹서

STANDARD SPECIFICATION

Model		PURIMIX
Size	Out(W x D x H)mm	132 x 185 x 130
	Weight(Kg)	2.3
Mixing	Motion	Orbital
Material		Powder Coated Aluminum
Speed Range		0 to 3,300 RPM
Speed Control		Stepless
Operating Mode		Touch or Continuous
Throw		Circular Orbit, 4.00mm diameter
Driving Motor		Brushless Shade Motor
Shaking Load Capacity		0.5Kg
Permissible Temp. & Humidity		5°C~40°C, < 80%
Power		AC 230V±10%, 50/60Hz
Power Consumption		30W
Application Container/Number		Hand holding; any micro tube or conical tube

FEATURE PLUS



3inch Rubber Platform Head

3인치 고무 플랫폼 헤드



Pop-Off Cup Head

Pop-off 컵 헤드



Remove Applicator

헤드 제거 툴

MAIN FEATURE



다기능 스위치

'연속사용 모드', '소프트 터치 모드', '전원오프'의 세 가지 모드를 쉽게 설정가능한 직관적인 스위치 차별화된 팝업 메커니즘을 구현한 '소프트 터치 모드'




강력한 퍼포먼스 및 안정성

강력한 모터와 견고한 금속 바디는 안정적인 믹싱을 가능케함



변환가능한 믹싱 속도

최대 3,300rpm 범위까지 믹싱 강도 조절이 가능함



다양한 어플리케이션

고무 헤드는 요철을 포함한 평면 디자인으로, 모양, 사이즈, 재질



높은 안정성

과열 및 과전류 방지 기능으로 모터를 보호함



오랜 수명의 모터

브러쉬가 없는 모터 장착



ORDER INFORMATION

구분	ITEM	ORDER CODE
		PURIMIX
제품	Main Body	PURIMIX-MB
스페이	3inch Rubber Platform Head	PMX-SP1
	Pop-Off Cup Head	PMX-SP2
	Remove Applicator	PMX-SP3